

616 Hollywood Avenue, Cherry Hill, NJ USA 1-800-333-3131


Werner-von-Siemens Str. 14a - 78224 Singen Germany +49 7731 912101

FINGER Instructions For Use

Become a fan of Keystone Industries

facebook.

visit us at www.keystoneind.com



INTENDED USE	USE WITHOUT GLOVE
The device consisting of the components described and shown in the illustration, was designed and manufactured exclusively to illuminate work areas close to hand (on which the device is fixed) that require artificial illumination. INAPPROPRIATE, REASONABLY FORSEEABLE USE In order to ensure safe use of the device any other use than the one described above is prohibited.	After having fixed the device it is possible to use it as follows (see illustration): 1) Switch on the LED 1-2 (Ref. 5) by pressing the "ON" button (Ref. 3); 2) Direct the luminous flux (Ref. 6) to the work area to be illuminated; 3) After finishing your work switch off the LEDs 1-2 (Ref. 5) by pressing the "OFF" button (Ref. 2); 4) Remove the device from your hand.
SAFETY WARNINGS	USE WITH GLOVE
 THE DEVICE IS NOT A TOY: KEEP IT OUT OF THE REACH OF CHILDREN; DO NOT USE THE DEVICE IN ENVIRONMENTS WITH DIFFERENT CHARACTERISTICS THEN THOSE LISTED IN THE TECHNICAL DATA; DO NOT SUBMERSE THE DEVICE IN LIQUIDS AND KEEP IT AWAY FROM SPLASHES DURING USE; DO NOT USE THE DEVICE UNLESS IT IS PROPERLY FIXED (AS SHOWN IN THE ILLUSTRATION); DO NOT USE THE DEVICE TO ILLUMINATE ROOMS; DO NOT USE THE DEVICE TO ILLUMINATE ROOMS; DO NOT ISE THE DEVICE TO ILLUMINATE ROOMS; DO NOT TOINT LED LAMPS DIRECTLY TOWARDS THE EYES; DO NOT TAMPER WITH ANY COMPORENT OF THE DEVICE: 	BE 1 HKM GMP : KGBE@ TO BE USED ONLY AND EXCLUSIVELY WITH LATEX OR NITRILE GLOVES CLASSIFIED AS MEDICAL EQUIPMENT (CE MARK), THE DEVICE IS TO BE CONSIDERED AS A CLASS I MEDICAL DEVICE AS PER REGULATION 12 OF DIRECTIVE 93/42/EEC. AS WELL AS MODIFIED DIRECTIVE 2007/47/EC. After having fixed the device it is possible to use it as follows (see illustration): 1) Switch on the LED 1-2 (Ref. 5) by pressing the "ON" button (Ref. 3); 2) Direct the luminous flux (Ref. 6) to the work area to be illuminated; 3) After finishing your work switch off the LEDs 1-2 (Ref. 5) by pressing the "OFF" button (Ref. 2)
DO NOT USE THE DEVICE WITH GLOVES THAT ARE NOT INTACT. ASSEMBLY (Battery replacement)	 4) Take off the glove and dispose of it in an appropriate container; 5) Remove the device from your hand.
Fix the device properly on the hand as follows (see illustration):	
FIXING THE BATTERY: 1) Open the battery compartment cover (Ref. 11); 2) Insert the battery supplied (Ref. 12) into the battery compartment (Ref. 10) observing the polarities as indicated; 3) Close the battery compartment (Ref. 11). FIXING THE CONTROL BOX: 1) Position the control box (Ref. 1) on your wrist; 2) Adjust the position of the wrist strap (Ref. 8) and fix it with the Velcro strip (Ref. 9);	WARRANTY The warranty is valid for 12 (twelve) months after the date of purchase and becomes void if the device is not used properly and/or not in accordance with the instructions and warnings described in the manual. The manufacturer does not assume any liability for damages resulting from failure, or partial operation resulting from any lack of conformity of the device. On signaling the lack of conformity the manufacturer shall communicate the possibilities for the collection and eventual replacement and/or repairs of the device, reserving the right to verify the cause of the malfunction.
FIXATION WITHOUT USING THE GLOVE:	DISPOSAL
 Insert the LED wires (Ref. 4) in the recesses at the inside of the two wire clamps (Ref. 7); Slip the two wire clamps (Ref. 7) onto your fingers (index and middle finger); Adjust the proper position of the LED wires (Ref. 4) so that the LEDs 1-2 (Ref. 5 a-b) are at the end of the finger. 	All components of the device must be identified according to the definitions of the EWC (European Waste Catalogue) and disposed of in strict compliance with the regulations of environmental protection currently in force.
FIXATION USING THE GLOVE:	KEYSTONE INDUSTRIES declares that the device
 FIXAIION USING I HE GLOVE: 1) Cut two pieces of adhesive tape (Ref. 15) of ideal length (circumstance of the fingers); 2) Adjust the proper position of the LED wires (Ref. 4) and fix them to the fingers (index and middle finger) using the two pieces of adhesive tape; 3) Put a new latex or nitrile glove on your hand in accordance with the use (Ref. 14). N.B. To replace the battery carry out steps 1) - 3). Dispose of the battery in accordance with the guidelines in chapter. 	AF 15 IOUE INUCUS INTES declares that the device denominated LED HAND ILLUMINATOR Mod. ET LIGHT Year of manufacture 2011 complies with directive EMC 2004/108/EC and to the following standards (il applicable) EN 55015 (2006) + A1 (2007) + A2 (2009) EN 61000-3-2 (2006) + A1 (2009) EN 61000-3-2 (2006) EN 61000-3-2 (2006)

INSTRUCTIONS (EN) - 8001235 / REV. 00



616 Hollywood Avenue, Cherry Hill, NJ USA *1-800-333-3131* Werner-von-Siemens Str. 14a - 78224 Singen Germany *+49 7731 912101*