



SAFETY DATA SHEET
C. E. T.® Enzymatic Toothpaste – Vanilla Mint Flavor

1. IDENTIFICATION

Product Name	C. E. T.® Enzymatic Toothpaste – Vanilla Mint Flavor
Recommended use of the chemical and restrictions on use	
Identified uses	Toothpaste for cats and dogs
Restrictions on Use	For veterinary use only
Product Numbers	CET201
Company Identification	Virbac AH, Inc. P.O. Box 162059 Fort Worth, Texas 76161
Customer Information Number	(800) 338-3659
Emergency Telephone Number	
CHEMTREC Number	(800) 424-9300
Other Emergency Number:	Poison Control Center: 1-800-222-1222
Issue Date	May 5, 2015
Supersedes Date	March 30, 2011

Safety Data Sheet prepared in accordance with OSHA's Hazard Communication Standard (29 CFR 1910.1200) and the Globally Harmonized System of Classification and Labelling of Chemicals (GHS)

2. HAZARD IDENTIFICATION

Hazard Classification
This product is classified as not hazardous in accordance with the Globally Harmonized System of Classification and Labelling (GHS).

Label Elements
Hazard Symbols
None

Signal Word: None

Hazard Statements
None

Precautionary Statements

Prevention
None

Response
None

Storage
None

Disposal
None

Other Hazards
None



2. HAZARD IDENTIFICATION

Specific Concentration Limits

The values listed below represent the percentages of ingredients of unknown toxicity.

Acute oral toxicity	<1%
Acute dermal toxicity	30 - 40%
Acute inhalation toxicity	50 - 60%
Acute aquatic toxicity	45 - 55%

3. COMPOSITION/INFORMATION ON INGREDIENTS

Synonyms:

This product is a mixture.

Component Name	CAS Number	Concentration
Sorbitol	50-70-4	30 - 40%
Amorphous silicon dioxide	112926-00-8	5 - 15%
Glycerine	56-81-5	1 - 10%
Titanium Dioxide	13463-67-7	0.1 - <1.0%

4. FIRST - AID MEASURES

Description of necessary first-aid measures

Eyes

Immediately flood the eye with plenty of water, holding the eye open. Obtain medical attention if soreness or redness persists.

Skin

If irritation develops wash skin thoroughly with soap and water. Obtain medical attention if redness or soreness persists.

Ingestion

Call a poison control center or doctor immediately for treatment advice. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person.

Inhalation

Remove person to fresh air. Seek medical attention if symptoms persist.

Most important symptoms/effects, acute and delayed

Aside from the information found under Description of necessary first aid measures (above) and Indication of immediate medical attention and special treatment needed, no additional symptoms and effects are anticipated.

Indication of immediate medical attention and special treatment needed

Notes to Physicians

Treat symptomatically.

5. FIRE - FIGHTING MEASURES

Extinguishing Media

Use extinguishing media appropriate for surrounding materials.



5. FIRE - FIGHTING MEASURES

Unusual Fire and Explosion Hazards

Can release hazardous vapors during a fire.

Protective Equipment for Fire-Fighting

Wear full protective clothing and self-contained breathing apparatus.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Wear appropriate protective clothing.

Environmental Precautions

Prevent the material from entering drains or watercourses.

Methods and materials for containment and cleaning up

Wipe up and transfer into suitable containers for recovery or disposal. Prevent the material from entering drains or watercourses.

7. HANDLING AND STORAGE

Precautions for safe handling

Wear appropriate protective clothing.

Conditions for safe storage

Store in original container in a cool, dry place. Store away from children and pets.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure limits are listed below, if they exist.

Glycerin (Mist)

ACGIH: TLV 10 mg/m³ 8h TWA.

OSHA: PEL 5 mg/m³ 8h TWA respirable fraction

15mg/m³ 8h TWA total dust

Silica: Amorphous, including diatomaceous earth

OSHA: PEL 20 mppcf 8h TWA

0.8 mg/m³ 8h TWA

The exposure limit is calculated from the equation, $80/(\%SiO_2)$, using a value of 100% SiO₂. Lower values of % SiO₂ will give higher exposure limits.

Titanium Dioxide

ACGIH TLV: 10 mg/m³ TWA

OSHA PEL: 15 mg/m³ TWA, total dust

Appropriate engineering controls

No specific measures necessary. Good general room ventilation is expected to be adequate to control airborne levels.



8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Individual protection measures

Respiratory Protection

Not required under normal conditions of use.

Skin Protection

Gloves

Eye/Face Protection

Not required under normal conditions of use.

Body Protection

Normal work wear.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance

Physical State	Solid (paste)
Color	White to off-white
Odor	Mint, sweet
Odor Threshold	No data available
pH	No data available
Density	No data available
Boiling Range/Point (°C/F)	No data available
Melting Point (°C/F)	No data available
Flash Point (PMCC) (°C/F)	Not flammable
Vapor Pressure	No data available
Evaporation Rate (BuAc=1)	No data available
Solubility in Water	No data available
Vapor Density (Air = 1)	No data available
VOC	No data available
Partition coefficient (n-octanol/water)	Not applicable
Viscosity	Not applicable
Auto-ignition Temperature	No data available
Decomposition Temperature	No data available
Upper explosive limit	No data available
Lower explosive limit	No data available
Flammability (solid, gas)	No data available

10. STABILITY AND REACTIVITY

Reactivity

Data is not available

Chemical Stability

Stable under normal conditions.

Possibility of hazardous reactions

Hazardous polymerization will not occur.

Conditions to Avoid

Heat - high temperatures



10. STABILITY AND REACTIVITY

Incompatible Materials

None known.

Hazardous Decomposition Products

Oxides of carbon - acrolein

11. TOXICOLOGICAL INFORMATION

Acute Toxicity

Sorbitol

Oral LD50 (rat) 15,900 mg/kg

Glycerin

Oral LD50 (rat) >5000 mg/kg

Dermal LD50 (guinea pig) >50,000 mg/kg

Inhalation LC50 (rat) >2.75 mg/L 4hr

Amorphous Silicon Dioxide

Oral LD50 (rat) > 5000 mg/kg

Dermal LD50 (rabbit) >2000 mg/kg

Specific Target Organ Toxicity (STOT) – single exposure

Sorbitol: Reports of adverse reactions to sorbitol are largely due to its action as an osmotic laxative when ingested orally, which may be exploited therapeutically. Ingestion of large quantities of sorbitol (> 20g/day in adults) should therefore be avoided.

Specific Target Organ Toxicity (STOT) – repeat exposure

Available data indicates this product is not expected to cause target organ effects after repeated exposure.

Serious Eye damage/Irritation

Available data indicates this product is not expected to cause eye irritation.

Skin Corrosion/Irritation

Available data indicates this product is not expected to cause skin irritation.

Respiratory or Skin Sensitization

Available data indicates this product is not expected to cause skin sensitization.

Carcinogenicity

Titanium Dioxide: IARC Overall Evaluation is 2B (Possibly carcinogenic to humans) IARC evaluation guidelines consider the generation of tumors, in 2 different studies within the same animal species, to be adequate criteria for an assessment of sufficient evidence. The conclusions of several epidemiology studies on more than 20000 TiO₂ industry workers in Europe and the USA did not suggest a carcinogenic effect of TiO₂ dust on the human lung. Mortality from other chronic diseases, including other respiratory diseases, was also not associated with exposure to TiO₂ dust. Based upon these studies, titanium dioxide is not expected to cause lung cancer or chronic respiratory diseases in humans at concentrations experienced in the workplace.

Germ Cell Mutagenicity

Available data indicates this product is not expected to be mutagenic.



11. TOXICOLOGICAL INFORMATION

Reproductive Toxicity

Available data indicates this product is not expected to cause reproductive toxicity or birth defects.

Aspiration Hazard

Not an aspiration hazard.

12. ECOLOGICAL INFORMATION

Ecotoxicity

Available data indicates this product is not expected to be ecotoxic.

Mobility in soil

No relevant studies identified.

Persistence/Degradability

No relevant studies identified.

Bioaccumulative Potential

No relevant studies identified.

Other adverse effects

No relevant studies identified.

13. DISPOSAL CONSIDERATIONS

Disposal Methods

Dispose of in accordance with all applicable local and national regulations.

14. TRANSPORT INFORMATION

Contact supplier for transport information.

15. REGULATORY INFORMATION

United States TSCA Inventory

This product contains ingredients that have not been verified for listing on the Toxic Substance Control Act Chemical Inventory.

Canada DSL Inventory

This product contains ingredients that have not been verified for listing on the Domestic Substance List (DSL).

California Proposition 65

This product contains the following materials which the State of California has found to cause cancer, birth defects or other reproductive harm: None

SARA Title III Sect. 311/312 Categorization

None



15. REGULATORY INFORMATION

SARA Title III Sect. 313

This product contains a chemical that is listed in Section 313 at or above de minimis concentrations. The following listed chemicals are present: None

16. OTHER INFORMATION

Legend

ACGIH: American Conference of Governmental Industrial Hygienists
BOD: Biological Oxygen Demand
CAS#: Chemical Abstracts Service Number
FIFRA: Federal Insecticide, Fungicide and Rodenticide Act
IARC: International Agency for Research on Cancer
LC50: Lethal Concentration 50%
LD50: Lethal Dose 50%
N/A: Denotes no applicable information found or available
NTP: National Toxicology Program
OSHA: Occupational Safety and Health Administration
PEL: Permissible Exposure Limit
STEL: Short Term Exposure Limit
TLV: Threshold Limit Value
TSCA: Toxic Substance Control Act

Revision Date: May 5, 2015
Replaces: March 30, 2011
Changes made: Update to GHS.

Information Source and References

This SDS is prepared by Hazard Communication Specialists based on information provided by internal company references.

Prepared By: EnviroNet LLC.

C.E.T. is a registered trademark of Virbac Corporation.

The information and recommendations presented in this SDS are based on sources believed to be accurate. Virbac AH, Inc. assumes no liability for the accuracy or completeness of this information. It is the user's responsibility to determine the suitability of the material for their particular purposes. In particular, we make NO WARRANTY OF MERCHANTABILITY OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, with respect to such information, and we assume no liability resulting from its use. Users should ensure that any use or disposal of the material is in accordance with applicable Federal, State, and local laws and regulations.