


# Covidien Needle and Syringe Selection Guide

## Monoject™ Veterinary Needles & Syringes



## Monoject™ Hospital Standard Needles & Syringes



## KenVet™ Veterinary Needles & Syringes



### The Choice is Yours

- Monoject™ Needles & Syringes, the industry leader
- KenVet™ Needles & Syringes, the economy line

### Choose Your Needle

- **Veterinary Angle Point:** Designed for strength to penetrate thick or dehydrated skin
- **Hospital Point:** Designed for smooth penetration into thin well hydrated skin
- **Metal insert** for robust needle attachment
- **Epoxy needle** attachment
- **Clear or colored** hub and sheath

### Choose Your Syringe

- **Luer lock or slip** tip syringe
- **Silver recycled or white** plunger

### Choose Your Packaging

- **Recyclable rigid** pack
- **Peel to open** soft pack

# Covidien Needle and Syringe Selection Guide

		Monoject™ Veterinary Needles & Syringes	Monoject™ Standard Needles & Syringes	KenVet™ Veterinary Needles & Syringes
Needle Sheath	Color Coded	■		
	Clear		■	■
Needle Point	Shallow needle angle point for strength to penetrate thick or dehydrated skin	■		
	Steep needle point angle for penetration into thin well hydrated skin		■	■
Needle Attachment Epoxy	Stabilizing metal insert	■		
	Hub attachment		■	■
Needle Hub	Clear	■		■
	Color coded		■	
Syringe Attachment	Luer lock	■	■	■
	Luer slip	■		
Syringe Graduations	Heat	■	■	
	Printed			■
Finger Flanges	Low profile	■	■	
	Extended			■
Plunger Color	Silver	■	■	
	White			■
Packaging	Rigid recyclable packaging	■	■	
	Soft pack packaging			■

[www.covidien.com/animalhealth](http://www.covidien.com/animalhealth)



COVIDIEN, COVIDIEN with logo, Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other Brands are trademarks of a Covidien company. ©2014 Covidien. All rights reserved.

H8202 xM 0214