


FELINE Flex-coil® Swivel Kit

Reorder No. FMCRE9003

U.S. Patent No. 5,163,923 • Canadian Patent No. 2,119,376

Sterile, disposable, I.V. system with 60 drop/ml drip chamber, roller clamp, swivel device, 3 foot polyurethane coil, on/off clamp, two Y-injection ports, releasable tie to attach set to patient, and 14 inch extension set with proximal female Luer lock connector and distal KINK-SAFE™ male Luer slip connector to reverse tubing direction after attaching to I.V. catheter.



INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171
e-mail: info@internationalwin.com • www.internationalwin.com