

# Foal/Calf I.V. Extension Set

## 24" Extension Set

Sterile, disposable, 24 inch extension set with small bore chew-resistant tubing, on/off clamp, proximal female Luer lock connector and distal KINK-SAFE™ male Luer slip connector to prevent kinking of tubing when the direction is reversed.



## INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171  
e-mail: [info@internationalwin.com](mailto:info@internationalwin.com) • [www.internationalwin.com](http://www.internationalwin.com)