


# Coiled Extension Set with Swivel

Reorder No. CE8020S

## Coiled Extension Set with Swivel

Sterile, disposable, extension set with 360 degree inline swivel, 10 foot polyurethane coil, 14 inches of tubing beyond coil, Y-injection site, on/off clamp for quick disconnect and restart, releasable strap to secure set to patient, proximal female Luer lock adapter and distal SPIN-LOCK™ male Luer connector. Coiled extension set usually used with infusion pumps attached to special primary I.V. sets.



## INTERNATIONAL WIN, LIMITED

340 North Mill Road, Suite 6, Kennett Square, PA 19348-2853 USA • (610) 444-0170 • (800) 359-4946 • Fax (610) 444-0171  
e-mail: [info@internationalwin.com](mailto:info@internationalwin.com) • [www.internationalwin.com](http://www.internationalwin.com)