



EnzClenz Enzyme Detergent

EnzClenz is a concentrated multi-enzyme detergent for cleaning surgical instruments, reusable diagnostic equipment and for general cleaning. All-purpose enzyme detergent is superior for cleaning and deodorizing kennels, work surfaces, walls & floors. Enzymes work to destroy animal blood, urine, feces and other organic material, rather than merely pushing waste around with soaps or masking odors with scents.

SKU#	DESCRIPTION	UOM	PACK
024176	Premise Disinfectants - EnzClenz Enzyme Detergent Concentrate	1 Gallon	BOTTL



Instrument Lubricant

Instrument Lubricant contains an all-natural corrosion inhibitor and lubricant specially formulated to lubricate and protect surgical instruments, hinges, box locks and general equipment during autoclaving.

Features: Biodegradable, phosphate free, chloride free, silicone free, rust inhibitor, compatible with all sterilization methods

SKU#	DESCRIPTION	UOM	PACK
014325	Instruments, Cleaning/Disinfecting - Lubricant Milk Bath	1 Gallon	BOTTL



Instrument Detergent

Used for cleaning glassware, instruments, plastic and rubber items.

Features: Biodegradable, phosphate free, chloride free, rust inhibiting, rinses spotless

SKU#	DESCRIPTION	UOM	PACK
015246	Instruments, Cleaning/Disinfecting - AP Pink Detergent	1 Gallon	BOTTL



MaxiCide with Activator

Used for the sterilization and high level disinfection of heat sensitive medical equipment.

Features: Glutaraldehyde base, contains a rust inhibitor, solution can be used for up to 28 days once mixed

SKU#	DESCRIPTION	UOM	PACK
038058	Instruments, Cleaning/Disinfecting - Maxicide with Activator 2.5% Buffered Glutaraldehyde	1 Gallon	BOTTL
038057	Instruments, Cleaning/Disinfecting - Maxicide Plus with Activator 3.4% Glutaraldehyde	1 Gallon	BOTTL



Maxizyme Tablets

Maxizyme cleaning tablets, one cleaner with multiple uses: ultrasonic cleaning, evacuation system cleaning, and instrument presoak. Each tablet contains protease enzymes that break down hard-to-clean deposits that chemical detergents leave behind. Gentle on instruments because of its nontoxic, neutral pH, biodegradable, noncorrosive formula that contains no chlorine or harsh alkaline. Easy to use, two tablets quickly dissolve in 1 gallon of warm tap water. Fresh, minty scent.

SKU#	DESCRIPTION	UOM	PACK
035890	Ultrasonic, Accessories - Maxizyme Tablets	64 Count	BOX



ULTRA Clean-L

Effectively cleans dried blood, tissue, pumice, food, rouge, Tripoli and buffing compounds from the surface of any appliance while providing a rust inhibitor.

SKU#	DESCRIPTION	UOM	PACK
037394	Ultrasonic, Accessories - Cleaner Solution ULTRA Clean-L	1 Gallon	BOTTL

PVL00103-1/9



Helping to improve your clinic's profitability, without lowering your quality of care.



Ultrasonic General Purpose Cleaner

A non-ammoniated concentrate. Effectively cleans dried blood, bone, tissue, oxide compounds, and more. Removes ions from instrument's surface to enhance the effectiveness of ultrasonic cleaning unit.

SKU#	DESCRIPTION	UOM	PACK
037646	Ultrasonic, Accessories - Cleaner Solution	1 Gallon	BOTTL

Autoclave Steam Indicator Tape

Autoclave Indicator Tape is designed to seal packs and provide visual evidence that packs have been exposed to the steam sterilization process.

Features: For sealing and indicating processed packages, adheres to all materials without leaving residue, for steam autoclave and chemical vapor



SKU#	DESCRIPTION	UOM	PACK
037281	Sterilization, Indicators - Indicator Tape Steam 0.5" x 60 yd	EACH	EACH
037282	Sterilization, Indicators - Indicator Tape Steam 0.75" x 60 yd	EACH	EACH
037283	Sterilization, Indicators - Indicator Tape Steam 1" x 60 yd	EACH	EACH

CSR Wrap

A protective wrap made from natural cellulose fibers for general-purpose sterilization. Ideally suited for operating rooms, central supply, kit packing and other areas where equipment is wrapped prior to sterilization. Material is strong and tear resistant, yet soft and flexible for easy wrapping. Allows excellent penetration for various methods of sterilization including steam, ethylene oxide, radiation, and E-beam. Autoclavable.



SKU#	DESCRIPTION	UOM	PACK
038029	Sterilization, Wraps - CSR Wrap 12" x 12"	100 Count	BOX
036045	Sterilization, Wraps - CSR Wrap 24" x 24"	100 Count	BOX
036046	Sterilization, Wraps - CSR Wrap 30" x 30"	100 Count	BOX
036039	Sterilization, Wraps - CSR Wrap 36" x 36"	100 Count	BOX

Self-Seal Dual Indicator Sterilization Pouch

Self-Seal Dual Indicator Sterilization Pouches with dual indicator offer superior quality and protection. These pouches feature the added convenience of internal and external indicators that are printed on the same side and change color when processed with steam and ethylene oxide.



SKU#	DESCRIPTION	UOM	PACK
037276	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 3.5" x 5.25"	200 Count	BOX
037274	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 3.5" x 9"	200 Count	BOX
037275	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 5.25" x 10"	200 Count	BOX
037277	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 7.5" x 13"	200 Count	BOX
037279	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 8" x 13"	100 Count	BOX
037278	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 10.5" x 17"	100 Count	BOX
037280	Sterilization, Pouches/Tubes - Autoclave Pouch Dual Indicator Self Seal 13" x 20"	100 Count	BOX



Helping to improve your clinic's profitability, without lowering your quality of care.