
### Insulin injection for pets

A STEP-BY-STEP GUIDE FOR PET OWNERS

### 

## Preparing the syringe

Be sure to read and follow the instructions that accompany the insulin prescribed by your veterinarian. For example, many types of insulin require gentle rolling of the vial to mix the contents before each use. Also make sure you're using the proper insulin syringe: U-40 and U-100 syringes are not interchangeable.



6.

5. If there are air bubbles inside the syringe: Gently tap on the side of the syringe, inject the bubbles back into the vial, and pull the plunger back to the prescribed dose.

Remove the needle from the vial. Now you're ready to give the injection to your pet. You need to keep the needle sterile, so don't let it touch anything.

# Giving the injection

Once the syringe has been filled with the correct dose of insulin, take a moment to calm your pet. Giving the injection will go more smoothly if both you and your pet are relaxed. With some pets, it may help to give the injection while they're eating. You may need to have someone help hold your pet so you can safely give the injection.

#### 1.

For dogs and cats, the most common places to inject insulin are on the neck and over the shoulder blades loose areas of skin are best. Follow your vet's recommendations and be sure to rotate injection sites so you aren't using the same location each time.

#### 2.

Gently pinch and pull up a section of your pet's skin to form a small "tent." Hold the syringe in your other hand and be sure not to let the needle touch anything.





#### 3.

Insert the needle through the skin into the base of the tent. Keep the needle approximately level with the surface of the skin so that it doesn't accidentally enter a muscle or go through the opposite side of the tent.

#### 4.

Place your thumb on the end of the plunger and gently push it all the way in to inject the entire dose of insulin.





#### 5.

Immediately pull the needle out of your pet's skin.


If you're not sure that your pet received the entire dose, do not inject more insulin. Wait until the next scheduled injection time.

It's a good idea to write down the date and time of each injection.

#### 6.

Safely store the used needle in an approved sharps container such as the UltiGuard<sup>®</sup> Safe Pack. Do not throw used needles into the trash.

Always follow state and local regulations when disposing of used needles.



Veterinarians can request more of these teaching aids by emailing us at reorder@ulticare.com. Please be sure to include your contact name, delivery address, and telephone number.

This information has been reviewed by Dr. Jane Armstrong, DVM, MS, MBA, Diplomate of the American College of Veterinary Internal Medicine. It is provided for educational purposes only and is not intended to be a substitute for professional and/or medical advice. For specific guidance on how to give an injection to your pet, talk to your veterinarian.



UltiCare® and UltiGuard® are registered trademarks of UltiMed, Inc. ©2014 UltiMed, Inc