



Veterinary Energy-based Systems





Veterinary Energy-based Systems

Covidien is the world leader in designing and manufacturing advanced energy-based medical systems and accessories worldwide. Our product offerings include the complete line of Covidien state-of-the-art electrosurgical generators and an extensive array of electrosurgical accessories. Whether your surgical requirements call for a cutting edge vessel sealing system for permanent vessel occlusion or for the perfect low power office generator with monopolar/bipolar capability, look no further than Covidien Animal Health.



ForceTriad™ Energy Platform

- Valleylab™ mode for a unique combination of monopolar hemostasis and dissection while using a lower power setting
- Instant Response™ technology for precise energy delivery to effectively seal the vessel or tissue bundle
- Next generation LigaSure™ system utilizes TissueFect™ technology for sealing vessels 2x faster than previous technology
- Bipolar instrument receptacle for a traditional bipolar instrument
- Monopolar output equipped with Return Electrode Monitoring (REM™) for reduced risk of patient burns

LigaSure™ Vessel Sealing System

- Instant Response™ technology for precise energy delivery to effectively seal the vessel or tissue bundle
- LigaSure™ tissue fusion technology for vessels and tissue bundles up to and including 7 mm
- Bipolar instrument receptacle for a traditional bipolar instrument




Force FX™ Electrosurgical Generator

- Instant Response™ technology for delivering consistent power as tissue impedance changes
- Monopolar output equipped with Return Electrode Monitoring (REM™) for reduced risk of patient burns
- Bipolar instrument receptacle for a traditional bipolar instrument



Tissue Response Technology vs Conventional Technology



Tissue Response Technology™ senses tissue impedance at the electrode site and adjusts current and output voltage 200 times per second to maintain a consistent surgical effect. This is the principle behind Covidien Instant Response™ technology.

- Tissue Response Technology
- Conventional Technology



SurgiStat™ II Electrosurgical Generator

- Bipolar instrument receptacle for a traditional bipolar instrument
- Monopolar output equipped with a return electrode contact quality monitoring system (RECCMS) for increased patient safety

Accessories



Force Triverse™ Pencil



REM Polyhesive™ II Patient Return Electrode



Ligasure™ Atlas



Buttonswitch Pencil



Bipolar Forceps

Order Information

Code	Description
FORCEFXC	Force FX™ FX-C Generator
FORCETRIAD	ForceTriad™ Energy Platform
SURGI20	SurgiStat™ II Generator
LIGASURE	LigaSure™ System Generator-Dom
FT3000	Force Triverse™ Electrosurgical Device 3m (10Ft)
E6008	Standard Footswitch
E2250H	Footswitch Pencil
E2350H	Rockerswitch Pencil
E2450H	Buttonswitch Pencil
E7507	REM Polyhesive™ II Patient Return Electrode

Code	Description
E751025	Infant Return Electrode
LS0300	LigaSure™ Single Footswitch
LS1037	LigaSure™ Atlas Laparoscopic SLR/DIV
LS1020	LigaSure™ Atlas 20cm
LS1200	LigaSure™ Precise
LS1500	LigaSure™ V Sealer/Driver
LS3090	LigaSure™ Max Instrument
LS 3091	LigaSure™ Max Disp Elect/Cord
LS3110	LigaSure™ TM XTD Instrument
LS3111	LigaSure™ TM XTD Disp Elec/Cord



For more information contact your local Covidien Animal Health sales representative or e-mail us at: animalhealth@covidien.com



COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or an affiliate.
©2008 Covidien AG or an affiliate. All rights reserved.

H6338 10M 09/08

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

1-800-962-9888
508-261-8000

WWW.COVIDIEN.COM