# Vedco, Inc 5503 Corporate Drive St. Joseph, MO 64507

SAFETY DATA SHEET

INQUIRY- (816) 238-8840 EMERGENCY- (800) 424-9300

Issue Date 03-Jun-2014 Revision Date 03-Jun-2014 Version 1

## 1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

Product identifier

Product Name D-256

Other means of identification

Product Code VED332 UN/ID No. 1903 Synonyms None

Recommended use of the chemical and restrictions on use

Recommended Use Veterinarian Grade Disinfectant Uses advised against No information available

Details of the supplier of the safety data sheet

**Manufacturer Address** 

Vedco, Inc. 5503 Corporate Drive, St. Joseph, MO 64507

Emergency telephone number

**Company Phone Number** 1-816-238-8840 **24 Hour Emergency Phone Number** 800-424-9300

Emergency Telephone Chemtrec 1-800-424-9300

#### 2. HAZARDS IDENTIFICATION

### Classification

### **OSHA Regulatory Status**

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Skin corrosion/irritation	Category 1 Sub-category B
Serious eye damage/eye irritation	Category 1

### Label elements

#### **Emergency Overview**

## Danger

## Hazard statements

Causes severe skin burns and eye damage

D-256



Physical state liquid **Odor** Floral Appearance Clear

#### **Precautionary Statements - Prevention**

Do not breathe dust/fume/gas/mist/vapors/spray Wash face, hands and any exposed skin thoroughly after handling Wear protective gloves and eye protection.

#### **Precautionary Statements - Response**

Immediately call a POISON CENTER or doctor/physician

Specific treatment (see .? on this label)

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing Immediately call a POISON CENTER or doctor/physician

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower

Wash contaminated clothing before reuse

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

Immediately call a POISON CENTER or doctor/physician

IF SWALLOWED: Rinse mouth. DO NOT induce vomiting

### **Precautionary Statements - Storage**

Store locked up

#### **Precautionary Statements - Disposal**

Dispose of contents/container to an approved waste disposal plant

#### Hazards not otherwise classified (HNOC)

## Other Information

• May be harmful if swallowed **Unknown Acute Toxicity** 

8.25% of the mixture consists of ingredient(s) of unknown toxicity

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No.	Weight-%	Trade Secret
Didecyldimethyl Ammonium Chloride	7173-51-5	<11	*

<sup>\*</sup>The exact percentage (concentration) of composition has been withheld as a trade secret.

#### 4. FIRST AID MEASURES

#### First aid measures

Eye contact Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids.

Consult a physician.

Skin Contact Rinse with clear water.

Inhalation Remove to fresh air. If breathing does not return to normal, seek medical attention.

Immediately drink large quantities of water. Get medical attention. Do not induce vomiting. Ingestion

Revision Date 03-Jun-2014

D-256

#### Most important symptoms and effects, both acute and delayed

**Symptoms** No information available.

#### Indication of any immediate medical attention and special treatment needed

Note to physicians Treat symptomatically.

### 5. FIRE-FIGHTING MEASURES

#### Suitable extinguishing media

Carbon dioxide (CO2). Dry chemical. Alcohol Foam.

Unsuitable extinguishing media Caution: Use of water spray when fighting fire may be inefficient.

#### Specific hazards arising from the chemical

No information available.

#### Explosion data

Sensitivity to Mechanical Impact None. Sensitivity to Static Discharge None.

#### Protective equipment and precautions for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

### 6. ACCIDENTAL RELEASE MEASURES

#### Personal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation, especially in confined areas.

Environmental precautions

**Environmental precautions** See Section 12 for additional ecological information.

#### Methods and material for containment and cleaning up

**Methods for containment** Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Dispose in accordance with federal and state regulations.

## 7. HANDLING AND STORAGE

#### Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice.

#### Conditions for safe storage, including any incompatibilities

**Storage Conditions** Keep out of reach of children. Keep container tightly closed in a dry and well-ventilated

place. Do not freeze. Follow label instructions.

Incompatible materials Strong acids. OXIDIZERS.

#### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Revision Date 03-Jun-2014

Floral

D-256

Control parameters

**Exposure Guidelines** This product, as supplied, does not contain any hazardous materials with occupational

exposure limits established by the region specific regulatory bodies.

**Appropriate engineering controls** 

**Engineering Controls** Showers

> Eyewash stations Ventilation systems.

#### Individual protection measures, such as personal protective equipment

Eye/face protection Tight sealing safety goggles.

Skin and body protection Wear Neoprene or protective rubber gloves. Eye wash station.

Respiratory protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved

> respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be

provided in accordance with current local regulations.

**General Hygiene Considerations** Do not eat, drink or smoke when using this product. Wash contaminated clothing before

### 9. PHYSICAL AND CHEMICAL PROPERTIES

#### Information on basic physical and chemical properties

Physical state liquid **Appearance** Clear Odor

Color **Odor threshold** No information available green

**Property** Values Remarks • Method

9.5 Hq

Melting point/freezing point

No information available 100 >212°F

Boiling point / boiling range

Flash point No information available

**Evaporation rate** 

Flammability (solid, gas)

Flammability Limit in Air

No information available

Upper flammability limit: No information available Lower flammability limit: No information available Vapor pressure No information available

Vapor density <1 **Specific Gravity** 1.0

Water solubility completely soluble Solubility in other solvents No information available Partition coefficient No information available **Autoignition temperature** No information available **Decomposition temperature** No information available Kinematic viscosity No information available Dynamic viscosity No information available **Explosive properties** No information available **Oxidizing properties** No information available

**Other Information** 

Softening point No information available Molecular weight No information available **VOC Content (%)** No information available

D-256 Revision Date 03-Jun-2014

Density 8.39

Bulk density No information available

## 10. STABILITY AND REACTIVITY

### Reactivity

No data available

### **Chemical stability**

Stable under recommended storage conditions.

#### **Possibility of Hazardous Reactions**

None under normal processing.

**Hazardous polymerization** Hazardous polymerization does not occur.

#### **Conditions to avoid**

Extremes of temperature and direct sunlight.

#### **Incompatible materials**

Strong acids. OXIDIZERS.

### **Hazardous Decomposition Products**

Carbon dioxide (CO2).

## 11. TOXICOLOGICAL INFORMATION

### Information on likely routes of exposure

Product Information No data available

**Inhalation** Irritation and difficulty in breathing.

**Eye contact** Severely irritating to eyes.

**Skin Contact** No data available.

**Ingestion** Gastric pain and vomiting.

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Didecyldimethyl Ammonium	= 84 mg/kg (Rat)	-	-
Chloride			
7173-51-5			

#### Information on toxicological effects

**Symptoms** No information available.

## Delayed and immediate effects as well as chronic effects from short and long-term exposure

Sensitization
Germ cell mutagenicity
Carcinogenicity
Reproductive toxicity
STOT - single exposure
STOT - repeated exposure
Aspiration hazard
No information available.

### Numerical measures of toxicity - Product Information

D-256 Revision Date 03-Jun-2014

**Unknown Acute Toxicity** 

8.25% of the mixture consists of ingredient(s) of unknown toxicity

The following values are calculated based on chapter 3.1 of the GHS document.

## 12. ECOLOGICAL INFORMATION

#### **Ecotoxicity**

27.53% of the mixture consists of components(s) of unknown hazards to the aquatic environment

#### Persistence and degradability

No information available.

### **Bioaccumulation**

No information available.

Other adverse effects No information available

## 13. DISPOSAL CONSIDERATIONS

#### Waste treatment methods

Disposal of wastes Disposal should be in accordance with applicable regional, national and local laws and

regulations.

**Contaminated packaging** Do not reuse container.

## 14. TRANSPORT INFORMATION

**DOT** Not regulated

## 15. REGULATORY INFORMATION

**International Inventories** 

TSCA Does not comply
DSL/NDSL Complies
EINECS/ELINCS Does not comply

ENCS Does not comply
IECSC Complies
KECL Does not comply
PICCS Complies
AICS Complies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances
IECSC - China Inventory of Existing Chemical Substances
KECL Kersan Existing and Evaluated Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances **PICCS** - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

#### **US Federal Regulations**

### **SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

### SARA 311/312 Hazard Categories

Acute health hazard	No
Chronic Health Hazard	No
Fire hazard	No
Sudden release of pressure hazard	No
Reactive Hazard	No

#### **CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

#### **CERCLA**

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material

#### **US State Regulations**

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals

## U.S. State Right-to-Know Regulations

## U.S. EPA Label Information

**EPA Pesticide Registration Number** Not Applicable

## **16. OTHER INFORMATION**

NFPA Health hazards 0 Flammability 0 Instability 0 Physical and Chemical

Properties -

Health hazards 2 Flammability 0 Physical hazards 0 Personal protection X

Prepared ByMoira BlackflowerIssue Date04-Jun-2014Revision Date10-Apr-2015

**Revision Note** 

No information available

### **Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.