



Jorgensen Laboratories, Inc.

A family owned business, owned by veterinarians for veterinarians, since 1965
(800) 525-5614, www.JorVet.com, Info@JorVet.com


Bone Binder

A great alternative – easy, intuitive, and fast...

Placing cerclage wire is a frustrating surgical procedure. They can break or loosen with time, plus are difficult to place securely. Cable ties are a simple and well known method of fixation. Their use in surgery has been limited because they have not been studied or approved as implants.

Braun/Aesculap has done extensive studying over time to demonstrate the safety of this product:

- Easy and fast application
- No loosening of Bone Binder
- No additional instruments
- Consists of histo compatible, non-absorbable Polyamide 6.6
- Contains very fine radiographic markers
- Each box contains 5 packs of 2 (10 total)
- Sterile



Three sizes:

- J1377A** Small. 2.5mm(W) x 90mm(L)
- J1377B** Medium. 3.5mm(W) x 120mm(L)
- J1377C** Large. 4.5mm(W) x 180mm(L)

