OLYMPIC FOLD-UP

Saves space 200 pound capacity 2 Lift-off stretchers



Ideal for clinics with limited space Moves patients in tight spaces and then folds up and stores out of the way.

An "extra" gurney and treatment table Also can be converted to a temporary exam-treatment table using the Hard Stretcher top.

Very lightweight and strong Weighs only 27 pounds, but can safely carry 200 pound dogs. Strong airplane construction with welded joints.

Unfolds and folds in seconds To unfold, just spread the sides apart and an automatic safety lock holds it open. To close for storage, release safety lock and push the sides together.

Rolls easily, turns quickly Two hand grips at both ends of gurney for easy steering in tight spaces. Large 5" casters roll smoothly. All casters lock.

2 Models Economy model of the Fold-Up comes with Soft Stretcher top. Standard model has 2 tops – a Soft and a Hard Stretcher.

Dimensions Footprint: Unfolded: 54"L x 18"W x 34"H. Folded: 54"L x 8"W.

3-year Warranty *All parts and labor,* except parts subject to normal wear, like fabric, casters, etc.





Compact storage Only 8-inches deep

Lift-off stretchers - safer for staff and patients





Interchangeable Olympic gurneys include 2 interchangeable stretcher tops. Stretchers can be easily lifted off, but are held securely by 4 posts.

Prevents back injuries Using the stretcher tops, even the heaviest patients can be transferred to tables, cages, or the floor without back injury.

Safer for patients Heavy patients are transferred on a stretcher, not hand carried. Less traumatic and more comfortable for patients.



Soft Stretcher For lifting and transferring patients. Plastic fabric is medical grade Herculite[®]. Anti-microbial. Wipes clean. X-ray transparent for gross evaluations.



Hard Stretcher Converts the Fold-Up into a treatment table. Molded of rigid reinforced polymer. Non-slip surface. Wipes clean. Drains into bucket (not furnished).

Ordering and Information

Olympic Fold-Up™. **Model 50604.** *Standard model.* With Soft Stretcher and Hard Stretcher. Fully assembled.

Olympic Fold-Up[™]. **Model 50602.** *Economy model.* With Soft Stretcher only. Fully assembled.