

# OLYMPIC FOLD-UP™

**Saves space 200 pound capacity 2 Lift-off stretchers**



**ROLLS**

**FOLDS**



**Compact storage**  
Only 8-inches deep

**Lift-off stretchers – safer for staff and patients**



Watch video demonstration  
[www.OlympicVeterinary.com](http://www.OlympicVeterinary.com) 

**Ideal for clinics with limited space** Moves patients in tight spaces and then folds up and stores out of the way.

**An “extra” gurney and treatment table** Also can be converted to a temporary exam-treatment table using the Hard Stretcher top.

**Very lightweight and strong** Weighs only 27 pounds, but can safely carry 200 pound dogs. Strong airplane construction with welded joints.

**Unfolds and folds in seconds** To unfold, just spread the sides apart and an automatic safety lock holds it open. To close for storage, release safety lock and push the sides together.

**Rolls easily, turns quickly** Two hand grips at both ends of gurney for easy steering in tight spaces. Large 5" casters roll smoothly. All casters lock.

**2 Models** Economy model of the Fold-Up comes with Soft Stretcher top. Standard model has 2 tops – a Soft and a Hard Stretcher.

**Dimensions** Footprint: Unfolded: 54"L x 18"W x 34"H.  
Folded: 54"L x 8"W.

**3-year Warranty** All parts and labor, except parts subject to normal wear, like fabric, casters, etc.

**Interchangeable** Olympic gurneys include 2 interchangeable stretcher tops. Stretchers can be easily lifted off, but are held securely by 4 posts.

**Prevents back injuries** Using the stretcher tops, even the heaviest patients can be transferred to tables, cages, or the floor without back injury.

**Safer for patients** Heavy patients are transferred on a stretcher, not hand carried. Less traumatic and more comfortable for patients.



**Soft Stretcher** For lifting and transferring patients. Plastic fabric is medical grade Herculite®. Anti-microbial. Wipes clean. X-ray transparent for gross evaluations.



**Hard Stretcher** Converts the Fold-Up into a treatment table. Molded of rigid reinforced polymer. Non-slip surface. Wipes clean. Drains into bucket (not furnished).

## Ordering and Information

**Olympic Fold-Up™. Model 50604. Standard model.** With Soft Stretcher and Hard Stretcher. Fully assembled.

**Olympic Fold-Up™. Model 50602. Economy model.** With Soft Stretcher only. Fully assembled.