

# OLYMPIC Hi-Lo WET TABLE™

**NEW**



- Always at a comfortable height
- Better patient presentation
- More efficient
- Faster procedures



**Watch video demonstration**   
[www.OlympicVeterinary.com](http://www.OlympicVeterinary.com)

**W**et tables are used every day in most clinics, but because they have no height adjustment, they are often uncomfortable and inefficient. Old-style wet tables cannot be adjusted to your height – whether you are taller or shorter – or for larger and smaller animals. One size wet table has to fit all.

The new Olympic Hi-Lo Wet Table has adjustable height – designed for operator comfort and for more efficient access to the patient during any type of procedure.

Just touch the foot switch and the electric lift raises and lowers the entire table. Now, you can work without back or shoulder strain, whether standing or seated. And whether the patient is a puppy or a large dog, the treatment site will always be at the right height for you – accessible and efficient.

This modern, adjustable wet table can save you valuable time. Procedures are faster and less tiring, so you can give more patients better care.

**Splash shield**  
Confines water.  
Prevents objects  
from falling  
behind table.



**Hand spray**  
6 ft. coil hose  
reaches length  
of table. Stores  
out of way.



**Extra large  
hair trap**  
Has a larger mesh  
surface, so it does  
not have to be  
cleaned as often.  
Easily lifts out of sink.

## Taller or shorter people can work without strain



**Table higher – 42"**  
A tall (6'1") person can raise the Hi-Lo Table to prevent back strain. No more extreme bending over to treat the patient. Reduces fatigue and provides better visualization of the treatment site.



**Table lower – 36"**  
People who aren't as tall (4'11") and are treating larger animals can adjust the height of the Hi-Lo Wet Table for the best patient access and a comfortable posture. More efficient, faster treatment.



**Seated with table raised – 34"**  
Ergonomic positioning for a taller person with a patient.



**Seated with table lowered – 32"**  
Comfortable for longer procedures and with all patients.

### Special Features

**Foot switch swivels** Up/down foot switch can turn to face you. Convenient and saves time. (See opposite page.)

**2 removable grates** Each stainless steel grate 22" x 22". Lift out for cleaning.

**Stainless steel exam top** Fits over grates to create a convenient exam surface. Faucet remains accessible when installed.

**Electrical outlets for accessories** 3 outlets on base of table.

**Toe room** Plenty of toe room on all sides of table when standing.

**Clean panels** Removable panels under base help keep hair and dirt out.

**Leveling feet** Adjustable to any floor.

**Dimensions** Table top: 60" x 24" including stainless tub, sink and faucets.

**Electrical** 120V/60Hz. All UL components. Table has ground fault safety breaker.

**Plumbing connections** Hidden at rear of table.

**Full 3-year Warranty** Includes all parts and labor.

### Can be installed in any clinic

The Olympic Hi-Lo Wet Table can easily be installed as a replacement for any older wet table or in a new location. It is designed to meet all local plumbing codes. Even though the table moves up and down, the exterior plumbing connections are stationary and can be mated to standard fittings for any plumbing systems. (For plumbing specifications, please inquire.)



Perpendicular to plumbed cabinet



Perpendicular to plumbed wall



Free standing when floor drain available

### Ordering and Information

**Olympic Hi-Lo Wet Table™. Model 51001.** Complete, ready to install. Includes faucet and removable grates. Designed to meet all local plumbing codes.

**Olympic Hi-Lo Wet Table™ with stainless steel exam top. Model 51002.** Includes faucet and removable grates. Designed to meet all local plumbing codes.