# Start Saving Today With The Pureline<sup>®</sup> M6000 Anesthesia Machine With Oxygen Concentrator



M6000

Stop paying for tanks, refills, and deliveries. Stop wasting time on tank billing, scheduling, management and maintenance. Instead, get your oxygen for free, the natural way—from the air





#### Savings for you, safety for your patients.

The Pureline M6000 anesthesia machine with an integrated oxygen concentrator is one of a kind! While similar machines have been in use in Europe for many years, Supera has perfected the technology, added safety features like the built-in E-tank manifold that allows mobile use of the machine without power. It also protects your patient in case of a power failure.

The Pureline oxygen concentrator is FDA approved, UL tested and backed by the industry's best warranty.

### Key Features

- Extracts FDA approved medical-grade oxygen (USP93%) from the air and delivers it on demand up to 5 lpm
- · Energy efficient
- · Built specifically for anesthesia machine use
- · Oxygen purity safety alarms
- Safety backup oxygen supply with integrated oxygen flush
- Made in the USA
- Accepts industry standard vaporizers
- Virtually maintenance free

P/N M6000 Includes safety pop-off valve P/N M6000NS Includes regular pop-off valve (not shown)



# The Pureline® M6000 Pays For Itself!

By eliminating the cost of oxygen tank rental, purchase, refills, delivery charges and tank management time, the Pureline M6000 will pay for itself very quickly.





Heavy-duty shelf holds all monitor types safely and securely.

# **Return On Investment**

can be in as little as 12 months!

Visit www.superavet.com/ROI to calculate YOUR actual savings!

#### The Benefits of the Pureline® M6000

- No tanks, hoses, regulators to buy, store or replace—ever
- (V)
- · Oxygen is "on demand"—never run out!
- · Accepts industry standard vaporizers
- 0-4 lpm precision flow meter
- Easy on/off; controls are large and easily accessible
- Stable and strong "A" frame chassis with a low center of gravity
- Small 18" square footprint saves valuable floor space
- Includes anesthesia accessories: adult and pediatric unilimb rebreathing circuits, non-rebreathing circuit and heavy-duty rebreathing bags
- Quiet 40 dBA operation (best in the industry)
- Includes one high pressure "E" tank for backup and oxygen flush (the back up e-tank is not filled by machine)

Learn more about how Pureline can benefit any veterinary practice on our website:

## www.superavet.com

Anesthesia Innovations Since 1993





Formerly LEI Medical (503) 723-5068 or (877) 620-1500 info@superavet.com 12330 S.E. Highway 212, Clackamas, OR 97015