Start Saving Today with the Pureline® OC6000 Anesthesia Oxygen Concentrator.



OC6000

Stop paying for tanks, refills, and deliveries. Stop wasting time on tank billing, scheduling, management, and maintenance. Instead, get your oxygen for free, the natural way—from the air.





Savings for you, safety for your patients.

Supera Anesthesia Innovations is proud to bring the revolutionary Pureline OC6000 oxygen concentrator to veterinarians in the U.S. for the first time. While similar machines have been in use in Europe for years, Supera has perfected the technology, added safety features, and now is making this money-saving equipment available here. The built-in E-tank manifold allows mobile use of the machine without power. It also protects your patient in case of a power failure. The Pureline OC6000 concentrator is FDA approved, UL tested, and backed by the industry's best warranty.

Key Features

- Extracts medical-grade oxygen (USP93%) from the air and delivers it on demand up to 5 lpm
- Connects with industry standard oxygen DISS fittings
- Safety backup oxygen supply
- Energy efficient; average 12 cents a day to operate
- Built specifically for anesthesia machine use
- Oxygen purity and flow safety alarms
- Made in the USA by craftsmen using state of the art manufacturing equipment
- Virtually maintenance free

The OC6000 pays for itself in just a few months

Want to see how the savings could pencil out for your practice? ▶





The Pureline® OC6000 Pays For Itself In months

By eliminating the cost of oxygen tank rental, purchase, refills, delivery charges and tank management time, the Pureline OC6000 will pay for itself very quickly.



can be in just months!

Visit www.superavet.com/ROI to calculate YOUR actual savings!





I.V. bag pole for attaching fluids, pumps or monitors



Heavy duty shelf for monitor or pumps. Built in organizers for bags or cords.



Industry standard DISS O2 connection. Auxiliary and Main supply switch is clearly marked and accessible.

The Benefits of the Pureline OC6000

- · No tanks, hoses, regulators to buy, store, or replace
- Oxygen is "on demand"—never run out!
- Quiet 40-dBA operation
- Stable and strong "A" frame chassis with a low center of gravity
- Small 18" square footprint saves valuable floor space
- Includes one high pressure "E" tank for backup and oxygen flush (the back up e-tank is not filled by machine)

Learn more about how Pureline can benefit any veterinary practice on our website:

www.superavet.com/pureline

Anesthesia Innovations Since 1993



Formerly LEI Medical

(503) 723-5068 or (877) 620-1500 toll-free

info@superavet.com

12330 S.E. Highway 212, Clackamas, OR 97015