
Presenting the OC6200 Oxygen Concentrator Convertible.



OC6200

Converts your <u>existing</u> anesthesia machine into a moneysaving oxygen concentrator anesthesia machine.





Key Features

Save money on tanks, refills, and delivery charges

- Save time on tank billing, scheduling, and management
- Extracts medical-grade oxygen (USP93%) from the air and delivers it on demand up to 5 lpm
- Energy efficient; average 12 cents/day to operate

concentrator—without replacing your existing equipment.

- Built specifically for anesthesia machine use
- Oxygen purity and flow safety alarms
- FDA approved, UL tested
- Industry's best warranty
- Sturdy, space-saving construction
- Made in the USA by craftsmen using state of the art manufacturing equipment
- · Virtually maintenance free

Adapters available for: Matrx VMS®, LEI Medical M1200, Shor-Line®, Supera M1200, CDS2000, CDS9000. Additional machine adapters are being developed, please call if your machine model is not listed.

Pureline® is a registered trademark owned by Lawson Enterprises, Inc. MATRX and VMS are registered trademarks owned by Midmark Corporation. All other trademarks are the property of their respective owners. Supera Anesthesia Innovations and the Pureline® OC6200 are not affiliated with nor endorsed by Midmark Corporation.



NOTE: Anesthesia machine and vaporizerine not included

Quality and Durability by Design.

The Pureline OC6200 works with your existing anesthesia machine. Your machine mounts securely on the durable Supera chassis and connects to the Pureline® oxygen concentrator. You start saving money immediately.

The Benefits of the Pureline M6000

- Premium-quality manufacturing
- Stable and strong "A" frame chassis with a low center of gravity
- Easy on/off; controls are large and easily accessible
- Oxygen purity and flow safety alarms
- Quiet 40-dBA operation
- Includes standard E-tank for backup—FREE!
- 3-in. heavy-duty locking casters
- Small 18" square footprint saves valuable floor space
- Large shelf with additional storage features
- Made in the USA by craftsmen using state of the art manufacturing equipment

Pays for itself in almost no time



Return On Investment

can be in just months!

Visit www.superavet.com/ROI to calculate YOUR actual savings!





Heavy Duty Shelf

- 17 ½" x 8" E
 - Built-in syringe holder
- Shown with M1200



Circuit and bag organizer Part # OC6220

- Holds four standard Unilimb circuits
- Three positions to hold rebreathing bags

Learn more about how Pureline can benefit any veterinary practice on our website:

www.superavet.com/pureline

Anesthesia Innovations Since 1993



Formerly **LEI Medical**

(503) 723-5068 or (877) 620-1500 toll-free

info@superavet.com

12330 S.E. Highway 212, Clackamas, OR 97015