## **Innovative Safety Products**



## Safety Pressure Relief Valve A revolutionary advance in anesthesia safety.

This patented pressure relief valve virtually eliminates the danger of over-inflating a patient's lungs.



With just a press of the pressure regulating button, you can instanly restrict the pressure to a safe level, while effectively ventilating the patient. *Safety for you and your patients* 





## **FEATURES:**

- Impossible to leave closed while unattended
- Normal pressure relief at 20-25 cmH20 using ventilation button
- Less than 1 cmH20 pressure at full open
- Less than 1cmH20 pressure at full closed
- Replaces any valve with 3/4-20 thread (the industry standard thread size)
- Proven life saving patented design
- 19mm hose connection
- 10-25 cmH20 pop-off pressure
- · Can be used with a ventilator

Safety Pressure Relief Valve P/N ACC311



www.SuperaVet.com

SUPERA

ANESTHESIA INNOVATIONS

(503) 723-5068 or (877) 620-1500 toll-free

info@superavet.com

12330 S.E. Highway 212 Clackamas, OR 97015