



Rescue Wipes One Step Disinfectant Cleaner & Deodorizer (US)

SECTION 1. IDENTIFICATION

Product Identifier Rescue Wipes One Step Disinfectant Cleaner & Deodorizer
Recommended Use Disinfectant Cleaner Wipes
Manufacturer Virox Technologies Inc., 2770 Coventry Rd., Oakville, ON, L6H 6R1, 905-813-0110
Emergency Phone No. Virox Technologies Inc., 1-800-387-7578
SDS No. 003543

SECTION 2. HAZARDS IDENTIFICATION

GHS Classification

Not classified under any GHS hazard classes.

GHS Label Elements

Signal Word: None.

Hazard Pictogram: None.

Hazard Statement(s): None.

Other Hazards

The product contains no substances which at their given concentration are considered to be hazardous to health.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Mixture:

Chemical Name	CAS No.	%	Other Identifiers
Hydrogen peroxide	7722-84-1	0.5	

Notes

Active ingredients are listed above. EPA Registration Number 74559-10

SECTION 4. FIRST-AID MEASURES

First-aid Measures

Inhalation

No specific first aid measures are required.

Skin Contact

No specific first aid measures are required.

Eye Contact

Flush with cool water. Remove contact lenses, if applicable, and continue washing.

Obtain medical attention if irritation develops or persists.

Ingestion

Not a normal route of exposure.

SECTION 5. FIRE-FIGHTING MEASURES

Extinguishing Media

Suitable Extinguishing Media

Not combustible. Use extinguishing agents compatible with product and suitable for surrounding fire.

Unsuitable Extinguishing Media

None known.

Specific Hazards Arising from the Chemical

None known.

Special Protective Equipment and Precautions for Fire-fighters

Wear self-contained breathing apparatus for fire fighting if necessary.

Product Identifier: Rescue Wipes One Step Disinfectant Cleaner & Deodorizer

SDS No.: 003543

Date of Preparation: February 08, 2016

SECTION 6. ACCIDENTAL RELEASE MEASURES

Personal Precautions, Protective Equipment, and Emergency Procedures

Use the personal protective equipment recommended in Section 8 of this safety data sheet.

Environmental Precautions

Before attempting clean-up, refer to hazard data. Review STORAGE and DISPOSAL section of label prior to disposal.

Methods and Materials for Containment and Cleaning Up

Never flush these wipes into toilets. Discard in solid waste bin. Towelettes contaminated with blood or body fluids should be disposed of according to federal, state and local regulations for infectious waste disposal.

SECTION 7. HANDLING AND STORAGE

Precautions for Safe Handling

Use good industrial hygiene practices in handling this material (see section 8). FOR COMMERCIAL AND INDUSTRIAL USE ONLY.

Conditions for Safe Storage

Store in an area that is out of direct sunlight. Avoid storage at elevated temperatures. KEEP OUT OF REACH OF CHILDREN.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Exposure Guideline

Chemical Name	ACGIH		OSHA PEL		AIHA WEEL	
	TWA	STEL	TWA	Ceiling	8-hr TWA	TWA
Hydrogen Peroxide	1 ppm		1 ppm			

Appropriate Engineering Controls

No specific ventilation requirements.

Individual Protection Measures

Eye/Face Protection

Not required if product is used as directed.

Skin Protection

Not required if product is used as directed.

Respiratory Protection

Not required if product is used as directed.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Basic Physical and Chemical Properties

Appearance	Clear colorless liquid saturated on wipes.
Odour Threshold	Not available
pH	2.0 - 2.5
Melting Point/Freezing Point	Not available
Initial Boiling Point/Range	Not available
Flash Point	> 200 °F (93 °C)
Evaporation Rate	Not available
Flammability (solid, gas)	Not applicable (liquid).
Upper/Lower Flammability or Explosive Limit	Not available (upper); Not available (lower)
Vapour Pressure	Not available
Vapour Density (air = 1)	Not available
Relative Density (water = 1)	1.006 at 20 °C
Partition Coefficient, n-Octanol/Water (Log Kow)	Not available
Auto-ignition Temperature	Not available
Viscosity	1.09 centistokes at 20 °C (kinematic)
Other Information	
Physical State	Wet wipes
Odour	Light Almond odour

Product Identifier: Rescue Wipes One Step Disinfectant Cleaner & Deodorizer

SDS No.: 003543

Date of Preparation: February 08, 2016

SECTION 10. STABILITY AND REACTIVITY

Reactivity

Not reactive.

Chemical Stability

This product is stable.

Possibility of Hazardous Reactions

Hazardous polymerization will not occur.

Conditions to Avoid

High temperatures.

Incompatible Materials

Do not mix with concentrated bleach products.

Hazardous Decomposition Products

None known.

SECTION 11. TOXICOLOGICAL INFORMATION

Likely Routes of Exposure

Inhalation; skin contact; eye contact; ingestion.

Acute Toxicity

LC50 (Inhalation): > 2.08 mg/L.

LD50 (Oral): > 5000 mg/kg.

LD50 (Dermal): > 5050 mg/kg.

Skin Corrosion/Irritation

Not classified under GHS criteria.

Serious Eye Damage/Irritation

Not classified under GHS criteria.

STOT (Specific Target Organ Toxicity) - Single Exposure

Inhalation

Not classified under GHS criteria.

Skin Absorption

Not classified under GHS criteria.

Ingestion

Not classified under GHS criteria.

STOT (Specific Target Organ Toxicity) - Repeated Exposure

Not classified under GHS criteria.

Respiratory and/or Skin Sensitization

Not classified under GHS criteria.

Carcinogenicity

Not classified under GHS criteria.

Reproductive Toxicity

Development of Offspring

Not classified under GHS criteria.

Sexual Function and Fertility

Not classified under GHS criteria.

Germ Cell Mutagenicity

Not classified under GHS criteria.

Interactive Effects

None known.

SECTION 12. ECOLOGICAL INFORMATION

This section is not required by OSHA.

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal Methods

Review the STORAGE and DISPOSAL instructions on product label prior to disposal.

Product Identifier: Rescue Wipes One Step Disinfectant Cleaner & Deodorizer

SDS No.: 003543

Date of Preparation: February 08, 2016

SECTION 14. TRANSPORT INFORMATION

Shipping Information:

Not regulated under Canadian TDG Regulations. Not regulated under US DOT Regulations.

Special Precautions Not applicable

for User

SECTION 15. REGULATORY INFORMATION

Safety, Health and Environmental Regulations

Canada

Domestic Substances List (DSL) / Non-Domestic Substances List (NDSL)

All ingredients are listed on the DSL/NDSL.

USA

Toxic Substances Control Act (TSCA) Section 8(b)

All ingredients are on the TSCA Inventory or are exempt from TSCA Inventory requirements under 40 CFR 720.

Additional USA Regulatory Lists

Other U.S. Federal Regulations

SARA 302/304/311/312 extremely hazardous substances: No listed substance.

SARA 302/304 emergency planning and notification: No listed substance.

US Regulations:

EPA Registration No.: 74559-10

This chemical is a pesticide product registered by the US Environmental Protection Agency and is subject to certain labelling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets (SDS) and for workplace labels for non-pesticide chemicals. The following is the hazard information as required on the pesticide label: KEEP OUT OF REACH OF CHILDREN.

California Proposition 65: This product is not subject to the reporting requirements under California's Proposition 65.

SECTION 16. OTHER INFORMATION

HMIS Rating Health - 0 Flammability - 0 Physical Hazard - 0

SDS Prepared By Virox Technologies Inc.

Phone No. (800) 387-7578

Date of Preparation February 08, 2016

Additional Information For an updated SDS please contact the supplier/ manufacturer listed on the first page of this document. Information contained herein was obtained from sources considered technically accurate and reliable. While every effort has been made to ensure full disclosure of product hazards, in some cases data is not available and is so stated. Since condition of actual product use are beyond control of the supplier, it is assumed that users of this material have been fully trained according to the requirement of all applicable legislation and regulatory instruments. No warranty, expressed or implied, is made and manufacturer/supplier will not be liable for any losses, injuries or consequential damages which may result from the use of or reliance on any information contained in this document. The contents of this document have been prepared in accordance with the OSHA Hazard Communication Standards.

Product Identifier: Rescue Wipes One Step Disinfectant Cleaner & Deodorizer

SDS No.: 003543

Page 04 of 04

Date of Preparation: February 08, 2016